

I may need to submit a review ?

What to do before submitting  
a case for review



**TRU-LINE™**  
INVISIBLE ORTHODONTIC SYSTEM

# 1. Check the latest aligner is fitting


Aligners need to fit snugly , if they do not treatment is not going to progress correctly. If the latest aligner is not fitting correctly try and back track a few of aligners and find one that fits correctly. Then have the patient wear this aligner for a week or 2 to try and catch up to following aligners. It is very useful to use the supplied Visual Track models to check to progress of the aligners.



VisualTRACK™

Visual Track Models will help identify the patients progress and compliance.

## 2. Check attachments



Attachments need to be clear from flash, this is vital for aligners to seat correctly. Also it's a good idea to polish the attachments but try not to distort the shape of the attachments, they are there to help the aligner seat correctly and to assist with special movements of particular teeth.

### 3. Loosen contacts



It will often be useful to loosen hard contact areas, and in fact it does not hurt to clean all contact areas. We suggest starting with a 0.1mm inter-proximal strip as this will help loosen, clean and reshape contact areas and allow teeth to move more freely. This may be needed more often on heavily crowded or difficult to move teeth. If teeth are in hard contact it will be difficult for them to move freely.

# 4. Aligner Chewies



Aligner Chewies help the aligners seat correctly and it is advised they be used for a few days when a new aligner is inserted. They can be invaluable when the aligner is not seating correctly, in fact in a day or so while using aligner chewies the aligner will fit into place correctly and tooth movement will commence

# 5. Check out the Truline Manual

In some circumstances other factors may effect tooth movement, or lack there of. A comprehensive list of general troubleshooting can be found in the TRULINE Manual. Contact us for a copy.

If you send a case for review you **MUST** to send us:

1. New impressions
2. Visual Track Models
3. Patient compliance record

\* Costs for revisions are clearly defined on the patients estimate. No complimentary revisions under Tru 10

**TRU-LINE**<sup>TM</sup>  
INVISIBLE ORTHODONTIC SYSTEM

**1300 878 336**

[www.trulinedental.com.au](http://www.trulinedental.com.au)

**FABDENT**<sup>TM</sup>  
Fabulous Products For Dentistry